




HDAC7A Antibody

Product Code	CSB-PA002916
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUI4
Immunogen	Synthesized peptide derived from the C-terminal region of Human Histone deacetylase 7a.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HDAC7; HDAC7A; Histone deacetylase 7; HD7; Histone deacetylase 7A; HD7a
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC7

Image



Western Blot analysis of HuvEc cells using Histone deacetylase 7a Polyclonal Antibody