




HKR1 Antibody

Product Code	CSB-PA002933
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10072
Immunogen	Synthesized peptide derived from the Internal region of Human HKR1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HKR1; ZNF875; Krueppel-related zinc finger protein 1; Protein HKR1; Zinc finger protein 875
Immunogen Species	Homo sapiens (Human)
Target Names	HKR1

Image



Western Blot analysis of COLO205 cells using HKR1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.