




HLA-DOA Antibody

Product Code	CSB-PA002934
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06340
Immunogen	Synthesized peptide derived from the Internal region of Human HLA-DO α .
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HLA-DOA; HLA-DNA; HLA-DZA; HLA class II histocompatibility antigen; DO alpha chain; MHC DN-alpha; MHC DZ alpha; MHC class II antigen DOA
Immunogen Species	Homo sapiens (Human)
Target Names	HLA-DOA

Image



Western Blot analysis of COLO205 cells using HLA-DO α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.