




HLA-DOA Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA002934 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P06340 |
| Immunogen | Synthesized peptide derived from the Internal region of Human HLA-DO α . |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | HLA-DOA; HLA-DNA; HLA-DZA; HLA class II histocompatibility antigen; DO alpha chain; MHC DN-alpha; MHC DZ alpha; MHC class II antigen DOA |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HLA-DOA |

Image



Western Blot analysis of COLO205 cells using HLA-DO α Polyclonal Antibody