




CHST10 Antibody

Product Code	CSB-PA002941
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43529
Immunogen	Synthesized peptide derived from the Internal region of Human HNK-1ST.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CHST10; Carbohydrate sulfotransferase 10; HNK-1 sulfotransferase; HNK-1ST; HNK1ST; HuHNK-1ST
Immunogen Species	Homo sapiens (Human)
Target Names	CHST10

Image



Western Blot analysis of HuvEc cells using HNK-1ST Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.