




HNRNPA2B1 Antibody

Product Code	CSB-PA002943
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22626
Immunogen	Synthesized peptide derived from the N-terminal region of Human hnRNP A2/B1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HNRNPA2B1; HNRPA2B1; Heterogeneous nuclear ribonucleoproteins A2/B1; hnRNP A2/B1
Immunogen Species	Homo sapiens (Human)
Target Names	HNRNPA2B1

Image



Western Blot analysis of HepG2 cells using hnRNP A2/B1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.