





HNRNPL Antibody

Product Code	CSB-PA002951
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14866
Immunogen	Synthesized peptide derived from the Internal region of Human hnRNP L.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HNRNPL; HNRPL; P/OKcl.14; Heterogeneous nuclear ribonucleoprotein L; hnRNP L
Immunogen Species	Homo sapiens (Human)
Target Names	HNRNPL

Image



Western Blot analysis of HuvEc cells using hnRNP L Polyclonal Antibody



Western Blot analysis of HELA cells using
hnRNP L Polyclonal Antibody