





SYNCRIP Antibody

Product Code	CSB-PA002953
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60506
Immunogen	Synthesized peptide derived from the Internal region of Human hnRNP Q.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SYNCRIP; HNRPQ; NSAP1; Heterogeneous nuclear ribonucleoprotein Q; hnRNP Q; Glycine- and tyrosine-rich RNA-binding protein; GRY-RBP; NS1-associated protein 1; Synaptotagmin-binding; cytoplasmic RNA-interacting protein
Immunogen Species	Homo sapiens (Human)
Target Names	SYNCRIP

Image



Western Blot analysis of SH-SY5Y cells using hnRNP Q Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using hnRNP Q Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.