





SYNCRIP Antibody

Product Code	CSB-PA002953
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60506
Immunogen	Synthesized peptide derived from the Internal region of Human hnRNP Q.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SYNCRIP; HNRPQ; NSAP1; Heterogeneous nuclear ribonucleoprotein Q; hnRNP Q; Glycine- and tyrosine-rich RNA-binding protein; GRY-RBP; NS1-associated protein 1; Synaptotagmin-binding; cytoplasmic RNA-interacting protein
Immunogen Species	Homo sapiens (Human)
Target Names	SYNCRIP

Image



Western Blot analysis of SH-SY5Y cells using hnRNP Q Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using
hnRNP Q Polyclonal Antibody