





HOXA11/HOXD11 Antibody

Product Code	CSB-PA002955
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31270/P31277
Immunogen	Synthesized peptide derived from the C-terminal region of Human HoxA11/D11.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HOXA11; HOX11; Homeobox protein Hox-A11; Homeobox protein Hox-11; HOXD11; HOX4F; Homeobox protein Hox-D11; Homeobox protein Hox-4F
Immunogen Species	Homo sapiens (Human)
Target Names	HOXA11/HOXD11

Image



Western Blot analysis of 293T-UV cells using HoxA11/D11 Polyclonal Antibody



Western Blot analysis of LOVO cells using HoxA11/D11 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.