




HOXD12 Antibody

Product Code	CSB-PA002964
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35452
Immunogen	Synthesized peptide derived from the C-terminal region of Human HoxD12.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HOXD12; HOX4H; Homeobox protein Hox-D12; Homeobox protein Hox-4H
Immunogen Species	Homo sapiens (Human)
Target Names	HOXD12

Image



Western Blot analysis of K562 cells using HoxD12 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.