




HGS Antibody

Product Code	CSB-PA002970
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14964
Immunogen	Synthesized peptide derived from Human Hrs around the non-phosphorylation site of Y334.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HGS; HRS; Hepatocyte growth factor-regulated tyrosine kinase substrate; Hrs; Protein pp110
Immunogen Species	Homo sapiens (Human)
Target Names	HGS

Image



Western Blot analysis of HuvEc cells using Hrs Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.