





HSPA8 Antibody

Product Code	CSB-PA002972
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11142
Immunogen	Synthesized peptide derived from the Internal region of Human HSC 70.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HSPA8; HSC70; HSP73; HSPA10; Heat shock cognate 71 kDa protein; Heat shock 70 kDa protein 8
Immunogen Species	Homo sapiens (Human)
Target Names	HSPA8

Image



Western Blot analysis of HeLa cells using HSC 70 Polyclonal Antibody



Western Blot analysis of Jurkat cells, NIH-3T3 cells, Rat muscle using HSC 70 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.