





HSPA8 Antibody

Product Code	CSB-PA002973
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11142
Immunogen	Synthesized peptide derived from the Internal region of Human HSC 70.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HSPA8; HSC70; HSP73; HSPA10; Heat shock cognate 71 kDa protein; Heat shock 70 kDa protein 8
Immunogen Species	Homo sapiens (Human)
Target Names	HSPA8

Image



Western Blot analysis of HeLa cells using HSC 70 Polyclonal Antibody



Western Blot analysis of A549 HELA cells using HSC 70 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.