




HSF1 Antibody

Product Code	CSB-PA002974
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00613
Immunogen	Synthesized peptide derived from Human HSF1 around the non-phosphorylation site of S303.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HSF1; HSTF1; Heat shock factor protein 1; HSF 1; Heat shock transcription factor 1; HSTF 1
Immunogen Species	Homo sapiens (Human)
Target Names	HSF1

Image



Western Blot analysis of COLO 293T cells using HSF1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.