




CDK5RAP3 Antibody

Product Code	CSB-PA002976
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JB5
Immunogen	Synthesized peptide derived from the Internal region of Human HSF-27.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDK5RAP3; IC53; MSTP016; OK/SW-cl.114; PP1553; CDK5 regulatory subunit-associated protein 3; CDK5 activator-binding protein C53; Protein HSF-27
Immunogen Species	Homo sapiens (Human)
Target Names	CDK5RAP3

Image



Western Blot analysis of 293 cells using HSF-27 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.