





DLGAP5 Antibody

Product Code	CSB-PA002993
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15398
Immunogen	Synthesized peptide derived from the C-terminal region of Human HURP.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DLGAP5; DLG7; KIAA0008; Disks large-associated protein 5; DAP-5; Discs large homolog 7; Disks large-associated protein DLG7; Hepatoma up-regulated protein; HURP
Immunogen Species	Homo sapiens (Human)
Target Names	DLGAP5

Image



Western Blot analysis of A549 cells using HURP Polyclonal Antibody



Western Blot analysis of K562 NIH-3T3 cells using HURP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.