





ICAM1 Antibody

Product Code	CSB-PA002997
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05362
Immunogen	Synthesized peptide derived from Human ICAM-1 around the non-phosphorylation site of Y512.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ICAM1; Intercellular adhesion molecule 1; ICAM-1; Major group rhinovirus receptor; CD antigen CD54
Immunogen Species	Homo sapiens (Human)
Target Names	ICAM1

Image



Western Blot analysis of HepG2 cells using ICAM-1 Polyclonal Antibody



Western Blot analysis of PC3 22RV1 HELA cells using ICAM-1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.