





IK Antibody

Product Code	CSB-PA003006
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13123
Immunogen	Synthesized peptide derived from the C-terminal region of Human IK.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IK; RED; RER; Protein Red; Cytokine IK; IK factor; Protein RER
Immunogen Species	Homo sapiens (Human)
Target Names	IK

Image



Western Blot analysis of HeLa cells using IK Polyclonal Antibody



Western blot analysis of HELA KB SH-SY5Y 293T 3T3 lysis using IK antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.