




IL10RA Antibody

Product Code	CSB-PA003017
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13651
Immunogen	Synthesized peptide derived from Human IL-10R α around the non-phosphorylation site of Y496.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL10RA; IL10R; Interleukin-10 receptor subunit alpha; IL-10 receptor subunit alpha; IL-10R subunit alpha; IL-10RA; CDw210a; Interleukin-10 receptor subunit 1; IL-10R subunit 1; IL-10R1; CD antigen CD210
Immunogen Species	Homo sapiens (Human)
Target Names	IL10RA

Image



Western Blot analysis of HepG2 cells using IL-10R α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.