





IL17RC Antibody

Product Code	CSB-PA003022
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NAC3
Immunogen	Synthesized peptide derived from the C-terminal region of Human IL-17RC.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL17RC; Interleukin-17 receptor C; IL-17 receptor C; IL-17RC; Interleukin-17 receptor homolog; IL17Rhom; Interleukin-17 receptor-like protein; IL-17RL; ZcytoR14
Immunogen Species	Homo sapiens (Human)
Target Names	IL17RC

Image



Western Blot analysis of HuvEc cells using IL-17RC Polyclonal Antibody



Western Blot analysis of LOVO cells using IL-17RC Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.