





IL1B Antibody

Product Code	CSB-PA003023
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01584
Immunogen	Synthesized peptide derived from the Internal region of Human IL-1 β .
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL1B; IL1F2; Interleukin-1 beta; IL-1 beta; Catabolin
Immunogen Species	Homo sapiens (Human)
Target Names	IL1B

Image



Western Blot analysis of HeLa cells using IL-1 β Polyclonal Antibody



Western Blot analysis of L929 MOUSE-BRAIN cells using IL-1 β Polyclonal Antibody

