




IL2 Antibody

Product Code	CSB-PA003024
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60568
Immunogen	Synthesized peptide derived from the Internal region of Human IL-2.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL2; Interleukin-2; IL-2; T-cell growth factor; TCGF; Aldesleukin
Immunogen Species	Homo sapiens (Human)
Target Names	IL2

Image



Western Blot analysis of HT29 cells using IL-2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.