




CSF2RB Antibody

Product Code	CSB-PA003029
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32927
Immunogen	Synthesized peptide derived from Human IL-3R β around the non-phosphorylation site of Y593.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSF2RB; IL3RB; IL5RB; Cytokine receptor common subunit beta; CDw131; GM-CSF/IL-3/IL-5 receptor common beta subunit; CD antigen CD131
Immunogen Species	Homo sapiens (Human)
Target Names	CSF2RB

Image



Western Blot analysis of HepG2 cells using IL-3R β Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.