





IL5 Antibody

Product Code	CSB-PA003033
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05113
Immunogen	Synthesized peptide derived from the Internal region of Human IL-5.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL5

Image



Western Blot analysis of HuvEc cells using IL-5 Polyclonal Antibody



Immunohistochemistry of paraffin-embedded human breast carcinoma tissue using CSB-PA003033 at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.