





ITGA3 Antibody

Product Code	CSB-PA003045
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26006
Immunogen	Synthesized peptide derived from the Internal region of Human Integrin α 3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ITGA3; MSK18; Integrin alpha-3; CD49 antigen-like family member C; FRP-2; Galactoprotein B3; GAPB3; VLA-3 subunit alpha; CD antigen CD49c
Immunogen Species	Homo sapiens (Human)
Target Names	ITGA3

Image



Western Blot analysis of HepG2 cells using Integrin α 3 Polyclonal Antibody



Western Blot analysis of 293 HepG2 RAT-STOMACH cells using Integrin $\alpha 3$ Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.