





ITGB3 Antibody

Product Code	CSB-PA003048
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05106
Immunogen	Synthesized peptide derived from Human Integrin β 3 around the non-phosphorylation site of Y785.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ITGB3; GP3A; Integrin beta-3; Platelet membrane glycoprotein IIIa; GPIIIa; CD antigen CD61
Immunogen Species	Homo sapiens (Human)

Image



Western Blot analysis of A549 cells using Integrin β 3 Polyclonal Antibody



Western Blot analysis of HepG2-UV MOUSE-BRAIN cells using Integrin β 3 Polyclonal Antibody