




IPMK Antibody

Product Code	CSB-PA003058
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFU5
Immunogen	Synthesized peptide derived from the C-terminal region of Human IPMK.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IPMK; IMPK; Inositol polyphosphate multikinase; Inositol 1; 3,4,6-tetrakisphosphate 5-kinase
Immunogen Species	Homo sapiens (Human)
Target Names	IPMK

Image



Western Blot analysis of COS7 cells using IPMK Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.