





IQGAP1 Antibody

Product Code	CSB-PA003059
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46940
Immunogen	Synthesized peptide derived from the N-terminal region of Human IQGAP1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IQGAP1; KIAA0051; Ras GTPase-activating-like protein IQGAP1; p195
Immunogen Species	Homo sapiens (Human)
Target Names	IQGAP1

Image



Western Blot analysis of HeLa cells using IQGAP1 Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using IQGAP1 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.