




# IQGAP3 Antibody

<b>Product Code</b>	CSB-PA003060
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q86VI3
<b>Immunogen</b>	Synthesized peptide derived from the N-terminal region of Human IQGAP3.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	IQ motif containing GTPase activating protein 3 antibody; IQGA3_HUMAN antibody; IQGAP 3 antibody; IQGAP3 antibody; MGC10170 antibody; MGC10831 antibody; MGC1947 antibody; OTTHUMP00000031854 antibody; OTTHUMP00000031855 antibody; Ras GTPase activating like protein IQGAP3 antibody; Ras GTPase-activating-like protein IQGAP3 antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	IQGAP3

## Image



Western Blot analysis of HuvEc cells using IQGAP3 Polyclonal Antibody

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.