





IRAK2 Antibody

Product Code	CSB-PA003062
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43187
Immunogen	Synthesized peptide derived from the Internal region of Human IRAK-2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IRAK2; Interleukin-1 receptor-associated kinase-like 2; IRAK-2
Immunogen Species	Homo sapiens (Human)
Target Names	IRAK2

Image



Western Blot analysis of 293 cells using IRAK-2 Polyclonal Antibody



Western Blot analysis of Jurkat HuvEc cells using IRAK-2 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.