




BAIAP2L1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA003069 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UHR4 |
| Immunogen | Synthesized peptide derived from the Internal region of Human IRTKS. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | BAIAP2L1; IRTKS; Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 1; BAI1-associated protein 2-like protein 1; Insulin receptor tyrosine kinase substrate |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | BAIAP2L1 |

Image



Western Blot analysis of HepG2 cells using IRTKS Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.