





TAOK3 Antibody

Product Code	CSB-PA003080
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2K8
Immunogen	Synthesized peptide derived from the Internal region of Human JIK.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TAOK3; DPK; JIK; KDS; MAP3K18; Serine/threonine-protein kinase TAO3; Cutaneous T-cell lymphoma-associated antigen HD-CL-09; CTCL-associated antigen HD-CL-09; Dendritic cell-derived protein kinase; JNK/SAPK-inhibitory kinase; Jun kinase-inhi
Immunogen Species	Homo sapiens (Human)
Target Names	TAOK3

Image



Western Blot analysis of Jurkat cells using JIK Polyclonal Antibody



Western Blot analysis of Jurkat cells using JIK Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.