




KCNH1 Antibody

Product Code	CSB-PA003093
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95259
Immunogen	Synthesized peptide derived from the Internal region of Human KCNH1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	KCNH1; EAG; EAG1; Potassium voltage-gated channel subfamily H member 1; Ether-a-go-go potassium channel 1; EAG channel 1; h-eag; hEAG1; Voltage-gated potassium channel subunit Kv10.1
Immunogen Species	Homo sapiens (Human)
Target Names	KCNH1

Image



Western Blot analysis of HuvEc cells using KCNH1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.