




ZNF18 Antibody

Product Code	CSB-PA003107
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17022
Immunogen	Synthesized peptide derived from the Internal region of Human KOX11.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF18; HDSG1; KOX11; ZKSCAN6; ZNF535; Zinc finger protein 18; Heart development-specific gene 1 protein; Zinc finger protein 535; Zinc finger protein KOX11; Zinc finger protein with KRAB and SCAN domains 6
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF18

Image



Western Blot analysis of HeLa cells using KOX11 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.