




ZNF24 Antibody

Product Code	CSB-PA003108
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17028
Immunogen	Synthesized peptide derived from the N-terminal region of Human KOX17.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF24; KOX17; ZNF191; ZSCAN3; Zinc finger protein 24; Retinoic acid suppression protein A; RSG-A; Zinc finger and SCAN domain-containing protein 3; Zinc finger protein 191; Zinc finger protein KOX17
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF24

Image



Western Blot analysis of HuvEc cells using KOX17 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.