





# KCNC4 Antibody

<b>Product Code</b>	CSB-PA003122
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q03721
<b>Immunogen</b>	Synthesized peptide derived from Human Kv3.4 around the non-phosphorylation site of S15.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Monkey
<b>Tested Applications</b>	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	KCNC4; Potassium voltage-gated channel subfamily C member 4; KSHIIIIC; Voltage-gated potassium channel subunit Kv3.4
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	KCNC4

## Image



Western Blot analysis of SH-SY5Y cells using Kv3.4 Polyclonal Antibody



Western Blot analysis of Jurkat cells using Kv3.4 Polyclonal Antibody

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.