




TOR1AIP1 Antibody

Product Code	CSB-PA003136
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JTV8
Immunogen	Synthesized peptide derived from the Internal region of Human LAP1B.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TOR1AIP1; Torsin-1A-interacting protein 1; Lamin-associated protein 1B; LAP1B
Immunogen Species	Homo sapiens (Human)
Target Names	TOR1AIP1

Image



Western Blot analysis of HeLa cells using LAP1B Polyclonal Antibody