




LILRA2 Antibody

Product Code	CSB-PA003155
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N149
Immunogen	Synthesized peptide derived from the Internal region of Human LIR-7.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LILRA2; ILT1; LIR7; Leukocyte immunoglobulin-like receptor subfamily A member 2; CD85 antigen-like family member H; Immunoglobulin-like transcript 1; ILT-1; Leukocyte immunoglobulin-like receptor 7; LIR-7; CD antigen CD85h
Immunogen Species	Homo sapiens (Human)
Target Names	LILRA2

Image



Western Blot analysis of HepG2 cells using LIR-7 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.