




AGPAT4 Antibody

Product Code	CSB-PA003161
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRZ5
Immunogen	Synthesized peptide derived from the Internal region of Human LPAAT- δ .
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AGPAT4; 1-acyl-sn-glycerol-3-phosphate acyltransferase delta; 1-acylglycerol-3-phosphate O-acyltransferase 4; 1-AGP acyltransferase 4; 1-AGPAT 4; Lysophosphatidic acid acyltransferase delta; LPAAT-delta
Immunogen Species	Homo sapiens (Human)
Target Names	AGPAT4

Image



Western Blot analysis of HeLa cells using LPAAT- δ Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.