





LRG1 Antibody

Product Code	CSB-PA003166
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02750
Immunogen	Synthesized peptide derived from the C-terminal region of Human LRG1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LRG1; LRG; Leucine-rich alpha-2-glycoprotein; LRG
Immunogen Species	Homo sapiens (Human)
Target Names	LRG1

Image



Western Blot analysis of hela cells using LRG1 Polyclonal Antibody



Western Blot analysis of Jurkat cells using LRG1 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.