




CACNG7 Antibody

Product Code	CSB-PA003170
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62955
Immunogen	Synthesized peptide derived from the C-terminal region of Human L-type Ca ⁺⁺ CP γ 7.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CACNG7; Voltage-dependent calcium channel gamma-7 subunit; Neuronal voltage-gated calcium channel gamma-7 subunit; Transmembrane AMPAR regulatory protein gamma-7; TARP gamma-7
Immunogen Species	Homo sapiens (Human)
Target Names	CACNG7

Image



Western Blot analysis of HepG2 cells using L-type Ca⁺⁺ CP γ 7 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.