





MAOA Antibody

Product Code	CSB-PA003189
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21397
Immunogen	Synthesized peptide derived from the Internal region of Human MAO-A.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAOA; Amine oxidase [flavin-containing] A; Monoamine oxidase type A; MAO-A
Immunogen Species	Homo sapiens (Human)
Target Names	MAOA

Image



Western Blot analysis of HeLa cells using MAO-A Polyclonal Antibody



Western Blot analysis of 22RV1 MCF7 cells using MAO-A Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.