




MARCKS Antibody

Product Code	CSB-PA003190
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29966
Immunogen	Synthesized peptide derived from Human MARCKS around the non-phosphorylation site of S163.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MARCKS; MACS; PRKCSL; Myristoylated alanine-rich C-kinase substrate; MARCKS; Protein kinase C substrate; 80 kDa protein, light chain; 80K-L protein; PKCSL
Immunogen Species	Homo sapiens (Human)
Target Names	MARCKS

Image



Western Blot analysis of HepG2 cells using MARCKS Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.