




SERPINB5 Antibody

Product Code	CSB-PA003195
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36952
Immunogen	Synthesized peptide derived from the Internal region of Human Maspin.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SERPINB5; PI5; Serpin B5; Maspin; Peptidase inhibitor 5; PI-5
Immunogen Species	Homo sapiens (Human)
Target Names	SERPINB5

Image



Western Blot analysis of HeLa cells using Maspin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.