




MC2R Antibody

Product Code	CSB-PA003204
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01718
Immunogen	Synthesized peptide derived from the C-terminal region of Human MC2-R.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MC2R; ACTHR; Adrenocorticotrop hormone receptor; ACTH receptor; ACTH-R; Adrenocorticotropin receptor; Melanocortin receptor 2; MC2-R
Immunogen Species	Homo sapiens (Human)
Target Names	MC2R

Image



Western Blot analysis of HeLa cells using MC2-R Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.