





MC5R Antibody

Product Code	CSB-PA003205
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33032
Immunogen	Synthesized peptide derived from the C-terminal region of Human MC5-R.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MC5R; Melanocortin receptor 5; MC5-R; MC-2
Immunogen Species	Homo sapiens (Human)
Target Names	MC5R

Image



Western Blot analysis of 3T3 cells using MC5-R Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using MC5-R Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.