





MERTK Antibody

Product Code	CSB-PA003228
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12866
Immunogen	Synthesized peptide derived from the N-terminal region of Human MerTK.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MERTK; MER; Tyrosine-protein kinase Mer; Proto-oncogene c-Mer; Receptor tyrosine kinase MerTK
Immunogen Species	Homo sapiens (Human)
Target Names	MERTK

Image



Western Blot analysis of HEPG2 cells using MerTK Polyclonal Antibody



Western Blot analysis of various cells using
MerTK Polyclonal Antibody