




MET Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA003231 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P08581 |
| Immunogen | Synthesized peptide derived from Human Met around the non-phosphorylation site of Y1234. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | MET; Hepatocyte growth factor receptor; HGF receptor; HGF/SF receptor; Proto-oncogene c-Met; Scatter factor receptor; SF receptor; Tyrosine-protein kinase Met |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MET |

Image



Western Blot analysis of HepG2 cells using Met Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.