





MGMT Antibody

Product Code	CSB-PA003238
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16455
Immunogen	Synthesized peptide derived from the N-terminal region of Human MGMT.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MGMT; Methylated-DNA--protein-cysteine methyltransferase; 6-O-methylguanine-DNA methyltransferase; MGMT; O-6-methylguanine-DNA-alkyltransferase
Immunogen Species	Homo sapiens (Human)
Target Names	MGMT

Image



Western Blot analysis of Jurkat cells using MGMT Polyclonal Antibody



Western Blot analysis of MCF7 PC3 cells using
MGMT Polyclonal Antibody